

Business Name: Dental Fix

Name: Mark Belmar

Address: 12340 Seal Beach Blvd, Suite B653

City: Seal Beach

State: CA

Zip: 90740

Phone: 7143307881

Fax:

Email: mark.belmar@dentalfixrx.com

Comments/Notes: According to the operation manual, a steady light indicates the system is overheating.

---

Equipment: Small Equipment Repair

---

Item #	Quantity	Make	Model	S/N	Problem
1	1	Dentsply	Cavitron Plus	G131-30052	"Service" light is on and nothing works.

---

Tracking Code: 1Z65VV921214928196