

Business Name: Doster Dental
Name: Beth Hearn RDH
Address: 1001 carter street ste b
City: Chattanooga
State: TN
Zip: 37402-0000
Phone: 4232658839
Fax:
Email: dentalcare@drdoster.net

Comments/Notes: We recently received this unit from you following repair, but it still does not operate correctly. Please repair and return ASAP. We appreciate your expedience as we need this unit operational for our RDH quickly. Thank you --Cross Doster DDS

Equipment: Handpiece Repair

Item #	Quantity	Make	Model	S/N	Problem
1	1	Parkell	TurboSensor Ultrasonic Scaler	54214	1. Too much water flow, 2. inconsistent response to foot pedal (power), 3. Missing front plate from power knob after repair

Tracking Code: 1Z65VV920316595005