

Business Name: Precision Fix, Inc.

Name: Mark Belmar

Address: 12340 Seal Beach Blvd, Suite B653

City: Seal Beach

State: CA

Zip: 90740

Phone: 7143307881

Fax:

Email: precision.fix1@gmail.com

Comments/Notes: Office reported water was leaking at the handpiece and where the handpiece cord goes into the unit. When I inspected the unit, I never saw water in those places. My guess is the water regulator is leaking.

Equipment: Small Equipment Repair

Item #	Quantity	Make	Model	S/N	Problem
1	1	Dentsply	Cavitron Plus	G131-30052	Water keeps dripping at tip.

Tracking Code: 1Z65VV921201515469